

WE ARE HERE TO HELP



We recognise that you may be facing new challenges within the healthcare system given the Coronavirus (COVID-19) outbreak. At Abbott, we have many ways to help you maintain contact with your patients using the FreeStyle Libre system during this difficult time.

FreeStyle Libre and Telehealth

If you are **unable to see your patients with diabetes** face to face in your clinic/practice, please remember you are able to access your patients' FreeStyle Libre glucose reports remotely via LibreView without them leaving their home or entering your clinic.

If you are **still seeing patients with diabetes in clinic**, we can support you in quickly getting them onto the FreeStyle Libre Digital Ecosystem so that if things do change, you will be able to access their glucose data on demand, remotely and securely. Our FreeStyle Libre Digital Ecosystem is free to access for the patient and healthcare professional.

NDSS Subsidy and FreeStyle Libre

The NDSS is still processing and approving eligibility forms for the new NDSS subsidy.

FreeStyle Libre does not require a starter kit which gives you more flexibility on how to manage patients starting through the NDSS subsidy.

Remote FreeStyle Libre education sessions: 1-on-1 remote education session with a patient or a remote group webinar for your patients

LibreView training: Training for yourself on how LibreView can assist you on monitoring your patients glucose management remotely

myINTERACT: Online education resources for use in remote consults. You may sign up for free access to myINTERACT through www.freestyleconnect.com.au

Group Starts: We will be in contact shortly regarding our FreeStyle Libre group start program and how we can assist you to continue delivering this.



**FreeStyle
Libre**
FLASH GLUCOSE MONITORING SYSTEM

Patients can do it remotely



life. to the fullest.

Abbott

LibreView, a product of Newyu, is supported on Windows and Mac-operating computers. Minimum system requirements for operating systems are Windows 7 or higher on Windows computers and Mac OS X Yosemite (10.10) or higher on Mac computers. Minimum system requirements for browsers are IE10 or higher, Chrome 37 or higher, Firefox 32 or higher, and Safari 8 or higher. It is recommended to use the latest version of your browser. LibreView can be viewed on iOS 6 or higher mobile devices. FreeStyle, Libre, and related brand marks are trademarks of Abbott Diabetes Care Inc. in various jurisdictions. Information contained herein is for distribution outside of the USA only. Abbott Diabetes Care, 666 Doncaster Road, Doncaster, Victoria 3108, Australia. If you have any questions, please contact our Customer Service Team 1800 801 478. ABN 95 000 180 389 © 2020 Abbott ADC-20398 V1.0

YOU CAN DO IT REMOTELY *LibreView*




FreeStyle
LibreLink

If you scan with the FreeStyle LibreLink app on your phone, you're already halfway there.

If you are using the FreeStyle LibreLink app, you **AUTOMATICALLY** have a LibreView account (you can access LibreView with the same log in details) and your data is **AUTOMATICALLY** uploaded to LibreView.



If you scan with the FreeStyle Libre reader, you need a LibreView account, **BEFORE** you can share your data.



- Go to www.LibreView.com
- Click on **Sign up**
- Choose the **LibreView account for patients** and follow on-screen instructions
- **Connect** your FreeStyle Libre reader
- Press **Upload data** and follow on-screen instructions


FreeStyle
Libre
FLASH GLUCOSE MONITORING SYSTEM



1800 801 478

If you have any questions or require support with the FreeStyle Libre system, please call FreeStyle Libre Customer Care or visit FreeStyleLibre.com.au



life. to the fullest.

Abbott

THIS PRODUCT MAY NOT BE RIGHT FOR YOU. READ THE WARNINGS BEFORE PURCHASE AND FOLLOW THE INSTRUCTIONS. Find out more at FreeStyleLibre.com.au

This product is indicated for measuring interstitial fluid glucose levels in people (age 4 and older) with insulin-requiring diabetes.

†Scanning the sensor to obtain glucose values does not require lancets. A finger prick test using a blood glucose meter is required during times of rapidly changing glucose levels when interstitial fluid glucose levels may not accurately reflect blood glucose levels or if hypoglycaemia or impending hypoglycaemia is reported by the System or when symptoms do not match the System readings.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC.

There are two easy ways to share your data with your healthcare professional

LibreView



Option A:

- 1 Ask your healthcare professional to **invite you** to share your data
- 2 You will receive an email inviting you to share. Click **Accept Invitation** and you will be directed to the **LibreView** website
- 3 Log into your **LibreView** account with the **same username and password** used to sign up for FreeStyle LibreLink.

If you don't use the FreeStyle LibreLink app and you don't have a LibreView account you will need to create an account first.

Option B:

- 1 Ask your healthcare professional for his/her **practice ID**
- 2 Log into LibreView or create a LibreView account. Then click on the **☰** icon, at the top right of the screen and select **Account Settings**
- 3 Click on **My Practices**
- 4 Enter your healthcare professional's Practice ID



1800 801 478

If you have any questions or require support with the FreeStyle Libre system, please call FreeStyle Libre Customer Care or visit FreeStyleLibre.com.au



life. to the fullest.

Abbott

THIS PRODUCT MAY NOT BE RIGHT FOR YOU. READ THE WARNINGS BEFORE PURCHASE AND FOLLOW THE INSTRUCTIONS. Find out more at FreeStyleLibre.com.au

This product is indicated for measuring interstitial fluid glucose levels in people (age 4 and older) with insulin-requiring diabetes.

†Scanning the sensor to obtain glucose values does not require lancets. A finger prick test using a blood glucose meter is required during times of rapidly changing glucose levels when interstitial fluid glucose levels may not accurately reflect blood glucose levels or if hypoglycaemia or impending hypoglycaemia is reported by the System or when symptoms do not match the System readings.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered

in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC.

11 The FreeStyle LibreLink app and the FreeStyle Libre reader have similar but not identical features. A finger prick test using a blood glucose meter is required during times of rapidly changing glucose levels when interstitial fluid glucose levels may not accurately reflect blood glucose levels or if hypoglycaemia or impending hypoglycaemia is reported by the FreeStyle LibreLink app. 12 The FreeStyle Libre sensor communicates with the FreeStyle Libre reader that started it or the FreeStyle LibreLink app that started it. A sensor started by the FreeStyle Libre reader will also communicate with the FreeStyle LibreLink app. 13 LibreLinkUp is a mobile application developed and provided by Newyu, Inc. Use of FreeStyle LibreLink

and LibreLinkUp requires registration with LibreView, a service of Abbott and Newyu Inc. FreeStyle Libre readers and sensors can be obtained from Abbott. 14 It is possible to accept the LibreLinkUp invitation and receive messages and warnings upon doing so or to reject the invitation. You should base your decision upon your knowledge and experiences in being able to react accordingly when obtaining a glucose value that is too high or too low. 15 Transfer of glucose data between applications depends on cellular connectivity. 16 LibreView data can be viewed in the Safari Browser on Mac OS X Mountain Lion or higher computers and on iOS 6 or higher mobile devices. Currently, uploading of glucose data is only supported on Windows-based computers. Minimum system requirements are Windows Vista with IE10 or the latest versions of Google Chrome and Mozilla Firefox, running on a 550MHz Pentium III, 512MB DRAM, 2G Hard Drive, USB 2.0, LCD Screen

with resolution of 1024X768. 17 LibreView is not intended to be a primary glucose monitor: home users must consult their primary device(s) and consult a healthcare professional before making any medical interpretation and therapy adjustments from the information in the software, and healthcare professionals should use information in the software in conjunction with other clinical information available to them.

FreeStyle, Libre, and related brand marks are trademarks of Abbott Diabetes Care Inc. in various jurisdictions. Information contained herein is for distribution outside of the USA only. Abbott Diabetes Care, 666 Doncaster Road, Doncaster, Victoria 3108, Australia. If you have any questions, please contact our Customer Service Team 1800 801 478. ABN 95 000 180 389 © 2020 Abbott ADC-20538 V1.0